

TRANSFORMATION OF *E. coli*

Before starting: pull LB pates with appropriate antibiotic out of the refrigerator to warm to room temperature and heat up water bath to 42°C.

1. Thaw DH_{5α} competent cells on ice for about 10 minutes.

(# of desired transformations + controls) x 100 μl = total μl of cells needed

2. In an ice bucket chill labeled 1.5 ml eppendorf tubes.
3. Aliquot 100 μl of thawed competent cells to chilled eppendorf tubes with wide tip pipette tips (blue tips)
4. Add 1 to 3 μl of plasmid DNA or ligation reaction to appropriate tube of cells and stir with tip gently. Keep tubes on ice for 30 minutes.
5. Heat shock the cells for 60 seconds in a water bath at exactly 42°C.
6. Gently transfer cells to ice for at least 1.5 minutes
7. Add 400 μl to 900 μl liquid LB to each tubes. Shake the tubes for 1 hour at 200 rpm and at 37°C. For Ampicillin resistant cells this 1 hour incubation step can be skipped.
8. Plate 50 μl and 200 μl of transformed bacteria by hockey-stick method. Plate negative controls first.
9. Incubate the plates in 37°C incubator for 14-16 hours.

Making plasmid glycerol stocks: Prior to plasmid DNA extraction, add equal amounts of 80% glycerol and bacterial cells in a 1.5 ml eppendorf tube and store in -80°C.