

Optimized by Veronica Vallejo and Ann Armenia  
March, 04

### **Preparing Competent Cells:**

#### **Before you start:**

500 ml LB broth	Autoclaved 250 ml cylinder
2 LB plates (NO AMP)	16 sterile orange cap tubes
2 sterile culture tubes	Orbital shaker at 37°C
Sterile spreaders	Incubator at 37°C
Sterile toothpicks	0.05 M CaCl <sub>2</sub> Autoclaved
Tweezers	80% glycerol autoclaved
Autoclaved 1 L flask	Autoclaved 10 ml glass pipetts

#### Day 1

1. From glycerol stock of DH5 $\alpha$  cells inoculate 2 LB plates. Note: We make two plates to ensure that we have at least one good colony to work with the next day.

#### Day 2

1. Place 5 ml of LB broth into two culture tubes. Inoculate single colony into each culture tube. Incubate overnight at 37°C shaking at 230 rpm. Note: Should grow between 16-20 hours
2. Autoclave 1 L flask and 250 ml cylinder and 500 ml LB media.

#### Day 3

1. Using the autoclave cylinder, place 250 ml of autoclaved LB broth into autoclaved 1 L flask. Place flask in incubator at 37°C to warm-up (~1 hr).
2. Add 5 ml overnight culture to the warmed 250 ml LB broth. Incubate at 37°C shaking at 230 rpm until A650 reads 0.25-0.3 (~2 h).
3. Place 10 50ml sterile orange cap tubes into -20°C to chill.
4. Place 0.05 M CaCl<sub>2</sub> into -20°C to chill.
5. Chill flask on ice and turn on the centrifuge (4°C).
6. Divide the culture into the 10 orange cap tubes (~25 ml of culture into each orange cap tube).
7. Centrifuge at 3500 rpm for 10 min at 4°C. Note: be sure to balance centrifuge. Place tubes on ice.
8. Decant supernatant, resuspend each pellet in 10 ml ice cold 0.05 M CaCl<sub>2</sub> by swirling gently.
9. Pool 2 cultures into 1 tube, such that you end up with 5 tubes total.
10. Place tubes on ice for 2 hours. Place ice bucket in 4°C to minimize melting.
11. Spin samples as before.
12. Decant supernatant, resuspend each pellet in a total volume of 5 ml or less of ice cold CaCl<sub>2</sub>, 15% glycerol. Note: for a total volume of 5 ml add: 937.5  $\mu$ l of 80% glycerol and 4.1 ml of ice cold CaCl<sub>2</sub>.
13. Place tubes on ice, then flash freeze in liquid Nitrogen and put immediately into -80°C freezer.